


UROLOGIE  
UROLOGY  
UROLOGIE  
UROLOGIA  
UROLOGIA



Empfohlene Operationsbestecke für den Bereich "UROLOGIE"  
(Details erhalten Sie über Ihren Fachhändler)

Suggested Instrument Sets for „UROLOGY“  
(Please ask your dealer for details)


Suggestions de boîtes opératoires pour le secteur "UROLOGIE".  
(Pour plus de détails veuillez contacter votre distributeur.)

Cajas quirurgicas para "UROLOGIA" (para detalles póngase en contacto  
con su distribuidor)


Strumentari raccomandati per "UROLOGIA"  
(Per ulteriori informazioni, vi preghiamo di contattare il Vostro rivenditore)

- SET 51.0100** Urology Basic Instrument Set
- SET 51.0200** Urology Paediatric Basic Instrument Set
- SET 51.0300** Bladder Operating Set (Basic)
- SET 51.1000** Urology – Stones Set
- SET 51.1500** Prostatectomy Set
- SET 51.1550** Perineal Prostatectomy Set
- SET 51.1600** Supra-Retropubic Prostatectomy Set
- SET 51.2000** Urethral Dilation Inter Urethrotomy Set
- SET 51.2100** Meatotomy Set
- SET 51.2200** Orchiectomy Set
- SET 51.2300** Cystoscopy Basic Set
- SET 51.2500** Epispadias Repair Set
- SET 51.3000** Vasectomy Set
- SET 51.3100** Vasectomy Set "Non-Scalpel"
- SET 51.5000** Circumcision Set - Adult
- SET 51.5100** Circumcision Set - Newborn


Prostata-Fasszange  
 Prostatic grasping forceps  
 Pince à saisir la prostate  
 Pinza para coger prostata  
 Pinza per afferare la prostata




**YOUNG**  
 51.0004  
 20 cm



**MILLIN**  
 23 cm




51.0009



51.0010


Kapselasszangen  
 Capsule forceps  
 Pincés pour saisir les capsules  
 Pinzas para coger capsulas  
 Pinza per afferare la prostata




**MILLIN**  
 51.0011  
 23 cm



für die vorderen Lappen  
 for anterior lobes  
 pour lobe antérieur  
 para lobo anterior  
 per lobuli anteriori



**MILLIN**  
 51.0015  
 23 cm



Fadenführungs-  
 zange  
 Ligatur carrier  
 Porte-ligature  
 Porta-ligaduras  
 Pinza porta filo

Blasenhalsspreizer  
Bladder neck spreader  
Ecarteur pour col vésical  
Dilatador para el cuello vesical  
Divaricatore per collo vescicale




**MILLIN**  
51.0020  
24 cm

Samenstrang-Fixierklemme  
Spermatic cord fixation clamp  
Pince pour fixer le canal déférent  
Pinza-fixajion canal espermatica  
Pinza per fissare cordone spermatico



**SPERMATIC**  
51.0040  
14,5 cm

Ureterklemmen  
Ureter clamps  
Pinces à urètere  
Pinzas clamps para ureter  
Pinze per uretere



**51.0042**  
21 cm #2

**LURZ**  
51.0041  
21 cm #1

Prostatahaken  
Prostatic Retractors  
Ecarteurs à prostate  
Sparadores para prostata  
Retractori prostatici




**LEGUEU**  
**51.0050**  
27 cm / 35 mm



**KOCHER**  
**51.0053**  
26 cm / 67 mm

Prostatahaken  
Prostatic Retractors  
Ecarteurs à prostate  
Sparadores para prostata  
Re Trattori prostatici



**YOUNG**  
**51.0055**  
21 cm

lateral  
lateral  
latéral  
lateral  
laterale



**YOUNG**  
**51.0056**  
22 cm

anterior  
anterior  
antérieur  
anterior  
anteriore



**YOUNG**  
**51.0057**  
21 cm

Bulbus-Haken  
Bulbus hook  
Crochet à bulbus  
Gancho para el bulbus  
Uncino per bulbo



**YOUNG**  
**51.0058**  
22 cm

gegabelt  
forked  
en fourhce  
horquilla  
a forchetta

Vasektomieklemmen / Vasectomy Clamps / Pincés pour Vasectomie / Pinzas para vasectomia / Pinze per Vasectomia



**51.0043**  
14 cm Ø 3,0 mm




**51.0044**  
14 cm Ø 4,7 mm




**51.0045**  
12,5 cm  
gebogen, spitz  
curved, sharp  
courbe, pointu  
curva, agudo  
curvo, acuto

Samenstranghaken / Spermatic cord hooks / Crochets pour cordon spermatique /  
Ganchos para cordón espermático / Uncini per cordone spermatico




**SPERMATIC**  
**51.0060**  
16 cm



**17.0001**  
16 cm



**MILLIN**  
**51.0020**  
28 cm



**MILLIN**  
**51.0080**

komplett, mit 2 Seitenvalven 57x25 mm und 80x25 mm, 1 Mittelvale 120x45 mm  
complete, with 2 lateral valves 57x25 mm and 80x25 mm, 1 center blade 120x45 mm  
complet, avec 2 valves latérales 57x25 mm et 80x25 mm, 1 valve centrale 120x45 mm  
completo, con 2 valvas laterales 57x25 mm y 80x25 mm, 1 valva central 120x45 mm  
completo, con 2 valve laterali 57x25 mm e 80x25 mm, 1 valva centrale 120x45 mm



**JUDD-MASSON**  
**51.0070**

Spreizweite 100 mm  
spread 100 mm  
ouverture 100 mm  
apertura 100 mm  
apertura 100 mm






## JUDD-MASSON

**51.0072**

komplett, mit 2 Seitenvalven 70x25 mm und 1 Mittelvale 84x24 mm  
complete, with 2 lateral valves 70x25 mm and 1 center blade 84x24 mm  
complet, avec 2 valves latérales 70x25 mm et 1 valve centrale 84x24 mm  
completo, con 2 valvas laterales 70x25 mm y 1 valva central 84x24 mm  
completo, con 2 valve laterali 70x25 mm e 1 valva centrale 84x24 mm



## HRYNTSCHAK

**51.0098**

Rahmen allein  
Frame only  
Cadre seul  
Bastidor solo  
Telaio solo

## HRYNTSCHAK

**51.0099**

Mittelvalve 35x100mm  
Center blade 35x100mm  
Valve centrale 35x100mm  
Valva central 35x100mm  
Valva centrale 35x100mm

## THOMSON-WALKER

**51.0090**

komplett, mit 2 Seitenvalven 70x25 mm und 75x25 mm, 1 Mittelvale 90x25 mm  
complete, with 2 lateral valves 70x25 mm and 75x25 mm, 1 center blade 90x25 mm  
complet, avec 2 valves latérales 70x25 mm et 75x25 mm, 1 valve centrale 90x25 mm  
completo, con 2 valvas laterales 70x25 mm y 75x25 mm, 1 valva central 90x25 mm  
completo, con 2 valve laterali 70x25 mm e 75x25 mm, 1 valva centrale 90x25 mm


## HRYNTSCHAK

**51.0097**

Spreizweite 110 mm  
spread 110 mm  
ouverture 110 mm  
apertura 110 mm  
apertura 110 mm

für Männer  
 for men  
 pour hommes  
 para hombres-  
 per uomini


fur Frauen  
 for women  
 pour femmes  
 para mujeres  
 per donne



- 51.0106 6 Fr.
- 51.0107 7 Fr.
- 51.0108 8 Fr.
- 51.0109 9 Fr.
- 51.0110 10 Fr.
- 51.0112 12 Fr.
- 51.0114 14 Fr.
- 51.0116 16 Fr.
- 51.0118 18 Fr.
- 51.0120 20 Fr.
- 51.0122 22 Fr.
- 51.0124 24 Fr.



- 51.0136 6 Fr.
- 51.0137 7 Fr.
- 51.0138 8 Fr.
- 51.0139 9 Fr.
- 51.0140 10 Fr.
- 51.0141 11 Fr.
- 51.0142 12 Fr.
- 51.0143 13 Fr.
- 51.0144 14 Fr.
- 51.0145 15 Fr.
- 51.0146 16 Fr.
- 51.0148 18 Fr.
- 51.0150 20 Fr.
- 51.0152 22 Fr.
- 51.0154 24 Fr.



51.0303  
 35 cm

**GUYON**  
**51.0302**  
 39 cm




Sekretlöffel biegsam  
 Secretion curette, malleable  
 Cuillère pour sécrétions, malleable  
 Cucharilla para secreción, maleable  
 Curetta per secrezioni, malleabile

**GUTTMAN**  
 51.0350  
 26 cm

Penisklemmen  
 Penis Clamps  
 Pincés pour pénis  
 Pinzas para pene  
 Pinze serrapene




**STOCKMANN**  
 51.0600




**STRAUSS**  
 51.0610

Beschneidungsinstrument  
 Circumcision instrument  
 Instrument de circoncision  
 Instrumento para circoncisión  
 Strumento per circoncisione




**MOGEN**  
 51.0720  
 7,5 cm



Fremdkörperzangen  
 Foreign body forceps  
 Pincés pour corps étranger  
 Pinzas para cuerpos extranjeros  
 Pinze per corpi estranei

**MATHIEU**  
 51.0800 20 cm  
 51.0802 28 cm




**VAN BUREN**

28 cm

- 51.1008** 8 Fr.
- 51.1010** 10 Fr.
- 51.1012** 12 Fr.
- 51.1014** 14 Fr.
- 51.1016** 16 Fr.
- 51.1018** 18 Fr.
- 51.1020** 20 Fr.
- 51.1022** 22 Fr.
- 51.1024** 24 Fr.
- 51.1026** 26 Fr.
- 51.1028** 28 Fr.
- 51.1030** 30 Fr.
- 51.1032** 32 Fr.
- 51.1034** 34 Fr.
- 51.1036** 36 Fr.

**51.1040** 8 - 36 Fr.

Satz / set / jeu / juego / serie



**GUYON**

26 cm

- 51.1107** 7Fr.
- 51.1108** 8Fr.
- 51.1109** 9Fr.
- 51.1110** 10 Fr.
- 51.1111** 11 Fr.
- 51.1112** 12 Fr.
- 51.1113** 13 Fr.
- 51.1114** 14 Fr.
- 51.1115** 15 Fr.
- 51.1116** 16 Fr.
- 51.1117** 17 Fr.
- 51.1118** 18 Fr.
- 51.1119** 19 Fr.
- 51.1120** 20 Fr.
- 51.1121** 21 Fr.
- 51.1122** 22 Fr.
- 51.1123** 23 Fr.
- 51.1124** 24 Fr.
- 51.1125** 25 Fr.
- 51.1126** 26 Fr.
- 51.1127** 27 Fr.
- 51.1128** 28 Fr.
- 51.1129** 29 Fr.
- 51.1130** 30 Fr.

**51.1140** 10 - 30 Fr.

Satz / set / jeu / juego / serie